

PATENT

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Group Art Unit: 1744

Examiner: Chin, Brad Y.

Atty. Dkt. No.: UTSC:669US

In re Application of:

Issam RAAD, Hend A. HANNA, and Nabeel

NABULSI

Serial No.: 10/044,842

Filed: January 11, 2002

For: NOVEL ANTISEPTIC DERIVATIVES

WITH BROAD SPECTRUM

ANTIMICROBIAL ACTIVITY FOR THE

IMPREGNATION OF SURFACES

DECLARATION OF ISSAM RAAD, HEND HANNA, AND NABEEL NABULSI UNDER 37 C.F.R. 8 1.131.

We, Issam Raad, Hend A. Hanna, and Nabeel Nabulsi, hereby declare as follows:

- We are the named inventors for the above-referenced patent application.
- 2. Prior to September 25, 1998, we conceived of the idea of the idea of preparing compositions that include a basic reagent and a dye, and methods for disinfecting or sterilizing a surface that involve applying to the surface a composition that includes a basic reagent and a dye.
- 3. As evidence of conception of the invention, attached as Exhibits 1-6 are copies a literature searches we conducted to assess what was known in the literature pertaining to certain anti-infective agents, two of which chlorhexidine and berberine. Our idea was to combine these

agents in a single composition, and to coat the surface of medical devices (such as central venous catheters) with these compositions in an effort to inhibit the growth of microbacterial organisms that cause device-related infections. The date of this search was prior to September 25, 1998. Berberine, one of the agents that we searched in our review of the literature, is a yellow plant dye. Chlorhexidine, another of the agents search in the literature review, is a basic reagent.

4. Furthermore, from prior September 25, 1998 until we filed our provisional application on January 12, 2001, we were diligent in conducting studies to prepare compositions of our invention and evaluate their effectiveness as antimicrobial compositions. During this period, there was continual activity on our lab on this project. As evidence of this activity, we provide Exhibits 7-16, which include additional literature searches for basic reagents and dyes (7-15), and a summary of experiments performed after September 25, 1998, but prior to January 12, 2001, which showed the efficacy of combining various basic reagents and dye (16). In this regard, references to "Gendine" are to a combination of Gentian violet (also noted at "Gv") and chlorhexidene), and references to "PCMX" are to [INSERT].

5. We hereby declare that all statements made by our own knowledge are true and all statements made on information and belief are believed to be true and further that statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment under § 100 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of this application or any patent issued thereon.

Date 11/14/05	TROSS
The state of the s	Issam Raad
Date	Hend A. Hanna
Date	
	Nabeel Nabulsi

Genidine [Gentian violet (GV) · Chlorhexidine (CHX)]

Impregnating with DCM solution of Genidine.

	MRSA ₂₀₆₆	Ps ₄₂₀₅	C. Parap. 1-100-
PVC (7.0 mm	28:28[24:26]	21:21[14:15]	0022 27:28[25:25]
Si ^{II}	17:17 (18:19) ^[] [21:20	5:3 (11:12) ^{[[} [6:0]	21:21(19:19) [23:24]
PU (2Lumen; 10 FR.)	21:21[19:19]	14:13*[14:14]	27:26[21:21]
Suture (silk)	17:17[16:16]	5:3[2:3]	21:21[12:12]

Immersed for 10 min.

	MRSA ₂₀₆₆	Ps ₄₂₀₅	C. Parap. ₁₋₁₀₀₋
PVC ^u (7.0 mm I.D)	28:29	22:23	27:27
Si (2 lumen, 10 FR)	19:19 (19:20) •	10:11 (12:13) ^{\$}	18:18 (24:25) ф
PU (2lumen; 10 FR.)	22:22	15:15	22:23
Suture	15:15	4:4	14:14

Immersed for 1h. Immersed for 2h.

Control 1 : CHX in (DCM + Methanol) and in Methanol

	MRSA	2066	Ps ₄	205	C. Parap.	C. Parap. 1-100-0022		
	DCM+MeOH	MeOH	DCM+MeOH	MeOH	DCM+MeOH	MeOH		
PVC	0:0	0:0	0:0	0:0	0:0	0:0		
Si	0:0	0:0	0:0	0:0	0:0	0:0		
PU	17:17	11:11	10:10	0:0	15:15	0:0		
Sutu	0:0	0:0	0:0	0:0	0:0	0:0		
re		l						

All immersed for 2 h.

Control 2: GV in DCM and methanol

I MDON		-	_
MRSA ₂₀₆₆	PS ₄₂₀₅	1 C.	Parap. ₁₋₁₀₀₋

Immersed for 2 h. except for the 10 FR. Silicone, which were immersed for 20 h.

Ovalues are for 5 FR. Single lumen and those in parenthesis are for a double-lumen Cook 10.0 FR catheter.

^{*}Gave a 17 mm zone against the multi-resistant PS4277, while minorifampin control yielded 3 mm.

values between [] are for addition of 2 eq. Base instead of 3 eq.

Cocond Twist.

 $^{^{\}phi}$ Values in parenthesis are for 20h immersion.

About 33% DCM/MeOH (v/v)

	<u> </u>				0.0	0022		
	DCM	MeOH [₱]	DCM	MeOH [₱]	DCM	MeOH [♠]		
PVC _{II}	25:25	20:21	0:0	0:0	27:27	18:19		
Si (2 lumen, 10 FR)	6:7	7:8	0:0	0:0	0:0	0:0		
PU (21umen; 10 FR.)	22:22	32:32	0:0	0:0	22:23	31:32		
Suture	8:8	10:11	0:0	0:0	0:0	0:0		

Immersed for 10 min. Immersed for 2 h.

^{\$\phi_{\text{All devices immersed for 2 h.}}\$}

Genidine in methanol:

	MRSA ₂₀₆₆	Ps ₄₂₀₅	C. Parap. ₁₋₁₀₀ -
PVC	24:25	13:13	23:23
Si	10:12	0:0	0:0
PU	17:17	7:0	16:17
Suture	10:10	0:0	5:6

Experimental:

A. Impregantion Procedure

7.35 ml of 1M solution potassium t-butoxide in THF was added to a solution of CHX diacetate (1.533g; 2.45 mmol) in 35 ml THF. The resulting heterogeneous solution was stirred for 20 min, then added to a solution GV (1.0 g; 2.45 mmol) in 30 ml THF. The mixture was stirred at ambient conditions for 1 h, then placed under the hood overnight to evaporate the solvent. The resulting residue was dissolved in 30 ml DCM. One-centimeter device segments were immersed in the DCM solution for the appropriate periods: PVC & PU for 10 min; Si & Silk Suture for 2 h. After removal of the devices from the solution, traces of solution was removed from the lumen, then placed under the hood to dry over night. The impregnated devices were washed with distilled water until the washings were colorless or very faint, then placed under an aseptic hood to dry under ambient conditions for at least 4 h, preferably over night.

B. Zones of Inhibition

BBL Mueller Hinton II agar plates were inoculated with 0.5 McFarland of the appropriate microorganism. The impregnated devices were embedded in the inoculated plates and placed in an incubator at about 37.5° for at least 18 h. Zones of inhibition were then measured and corrected for yeast after incubation for several additional hours.

Durability of Polyurethane & Silicone Impregnated with GV*CHX

Zones of Inhibition against MRSA₂₀₆₆.

Tones of	r rimirn	TCTOH S	igarnsc	PIKSH206	6 •			
Catheter	Day=0	Day=3	Day=7	Day=10	Day=17	Day=30	Day=45	Day=61
PU	21:21	18:18	16:16	15:14	14:14	11:11	8:8	
Si	21:21	18:18	15:15	11:11	11:11	7:8	0:0	
	Day=75	********			alex needs to see the			
PU ·								
Si								

Alternative preparation of Gendine

The neutral form of chlorhexidine [55-56-1] is used instead of the salt form. Hence, .0025 mol of chlorhexidine is added to a stirring heterogeneous solution of 1 g of GV in 60 ml anhydrous THF at room temperature, and the resulting mixture is stirred for 1 h, then placed under the hood to evaporate the solvent. The resulting residue was dissolved in 30 ml DCM. The product did not totally dissolve in DCM.

	MRSA ₂₀₆₆		Ps	4205	C. Parar	C. Parap. 1-100-0022		
	Salt	Neutral	Salt	Neutral	Salt	Neutral		
PVC	28:28	29:29	21:21	19:19	27:28	30:30		
Si	18:19	21:21	11:12	9:9	21:21	22:22		
PU	21:21	23:23	14:13	12:12	27:26	22:22		
Sutu re	17:17	16:16	5:3	2:4	21:21	15:15		

Genidine in methanol:

Again placed under the hood to evaporate DCM. The resulting residue was dissolved in 30 ml MeOH, but the residue did not totally dissolve.

	MRS	A ₂₀₆₆	Ps	4205	C. Parap. ₁₋₁₀₀₋		
	Salt	Neutral	Salt	Neutral	Salt	Neutral	
PVC	24:25	21:22	13:13	12:12	23:23	20:20	
Si	10:12	15:15	0:0	0:0	0:0	8:9	
PU	17:17	19:19	7:0	12:12	16:17	12:12	
Suture	10:10	12:13	0:0	0:0	5:6	7:7	

Durability of Gendine-Coated Silicone UT Catheter

Zones of Inhibition against EN3836(VRE and E. Coli 3226 after immersion in Urine.

Organism	Day=0	Day=3	Day=7	Day=14	Day=21	Day=28	Day=35	Day=42	Day=49
EN3836	23:23	20	19:19	17:15	13:15	14:13	13		
E. Coli	18:18	15	13:13	13:12	12:12	12:11	11		

Day=0 against MRSA₂₀₅₆ = 26:27; against Ps = 18:18; against C. parap. = 24:25

Durability of Gendine-Coated ET PVC Tube

Zones of Inhibition against PS after immersion in Urine.

EN ₃₈₃₆ 23:23 20 19:19 17:15 13:15 14:13 13	Organism	Day=0	Day=3	Day=7	Day=14	Day=21	Day=28	Day=35	Day=42	Day=49
	EN3836	23:23	20	19:19	17:15			13		
				0.000						

Impregnating with Gendine in BuOAc at 40° C

The procedure is that adopted from the patent. CHX was added to a solution of GV in 30 ml n-BuOAc at room temp, then about 3 ml of MeOH was added. The resulting mixture along with the devices were heated for one hour at 40°C.

	MRSA ₂₀₆₆	Ps4205	C. Parap. ₁₋₁₀₀₋₀₀₂₂
PVC	27:27	17:17	26:26

Si	17:17	0:0	18:18	
PU	19:19	13:14	20:20	
Silk	14:14	5:5	14:14	

Impregnating with Gendine in 25 ml DCM + 5 ml BzOH

	MRSA ₂₀₆₆	Ps4205	C. Parap. ₁₋₁₀₀₋₀₀₂₂
PVC	23:23	13:17	23:23
Si	18:18	0:0	17:17
PU	20:20	17:17	20:20
Silk	11:11	0:0	9:9

Impregnating Biliary Stents with Gendine

One-centimeter segments were immersed in a DCM solution of Gendine over night. After drying overnight at ambient conditions, the segments were placed in a tube, washed with distilled water, then dried under the aseptic hood over night at ambient conditions. Then the impregnated pieces were embedded in inoculated Mueller-Hinton agar plates, and incubated over night. The resulting zones are given below.

Zone of Inhibition (mm) for Gendine-

impregnated Biliary Stents.

Microorganism	Zone
MRSA ₂₀₆₆	20
Ps4205	8
E. Coli ₃₂₀₂	10:11
E. Coli ₃₂₀₃	11:11
E. Coli ₃₂₂₆	11
Kb ₂₄₆₁	9:9
Kb ₂₅₄₈	9:10
Kb ₂₅₅₆	6:10
EN3836(VRE)	17
C. Albican ₆₄₅₅₁	25
C. Parap ₁₋₁₀₀₋₀₀₂₂	18

The molecular structures and electronic states of gentian violet has been the subject of extensive studies driven by the observed inhomogeneity of its absorption spectra. In agreement with our observations, Goldacre and Phillips demonstrated the bleaching of gentian violet in the presence of hydroxide. They attributed bleaching to nucleophilic attack of the hydroxide unto the benzylic carbonium center. 2

Similar bleaching is also observed in our laboratory for Gendine, both in methanol and dichloromethane. This is consistent with the presence of both uncharged and a charged gentian violet moiety in Gendine, where the latter is responsible for imparting the color

due to the extensive electronic delocalization. Consequently, and similar to that found for the acridinium dye Acriflavin³, Gendine must exists initially as an EDA (electron Donor-

$$\begin{bmatrix} \delta^+ & \delta^+ & \delta^+ & \delta^+ \\ \delta^+ & \delta^+ & \delta^+ & \delta^+ \\ \delta^+ & \delta^+ \\ \delta^+ & \delta^+ \\ \delta^+ & \delta^+ & \delta^+ \\ \delta^+ & \delta^$$

Acrilflavin

Acceptor) complex formed between the cationic gentian violet and chlorhexidine, which is also similar to that observed for gentian violet and pyridine (GV⁺·Py). ⁴

As a result, and at any given time, a solution of Gendine can consist of a mixture of both structural isomers. Meanwhile, zones of inhibition exhibited by devices impregnated with

Dark blue solution

Golden Solution

both solutions rule out conformational isomerism and confirm the presence of structural isomers, since it is expected for conformational isomers to impart similar zones. The data given below support existence of the two structural isomers for Gendine, the EDA "charge-transfer" complex and the covalently bonded isomer.

Zones of Inhibition (mm) for Gendine-Impregnated Devices

-	Blue Solution			Golden Solution			
	MRSA	PS	C. Parap.	MRSA	PS	C. Parap.	
PVC	28:28	22:22	27:28	22:23	15:15	21:21	
PU	21:21	15:15	22:23	19:19	13:13	17:17	

¹The color of the device is dark blue when impregnated with the blue solution, and is light gold (turns darker over night) when impregnated with the light golden solution.

II. Durability and Stability

For long-term intravenous therapy, Schierholz et al. pointed out the importance of continued release (exceeding 10 days) of antimicrobial, and attributed failure of infection prevention with chlorhexidine- silver sulfadiazine coated catheters to the decreased release beyond 48 hours.⁵

For gendine-impregnated catheters, the efficacy and stability in human serum is being studied for polyurethane (PU) and silicone (Si). As of this date, results show continued release of gendine bevond 30 days (Table 1).

Table 1. Zone of Inhibition (mm) against MRSA₂₀₆₆ after incubation in Human Serum

Catheter	Day=0	Day=3	Day=7	Day=10	Day=17	Day=31	Day=45	Day=60
PU	21:21	18:18	16:16	15:14	14:14	- 11:11		
Si	21:21	18:18	15:15	11:11	11:11	7:8		

References

The color of the device is dark blue when impregnated with the blue solution, and is light gray when impregnated with the light golden solution.

¹ Maruyama, Y.; Ishikawa, M.; Satozono, H. J. Am. Chem. Soc. 1996, 118, 6257-6263.

² Goldacre, R. J.; Philips, J. N. J. Chem. Soc., 1949, 1724-32.

³ Merck Index 12, 125.

⁴ Liang, E. J.; Ye, X. L.; Kiefer, W. J. Phys. Chem., 1997, 101, 7330-7335.

⁵ Schierholz, J.; Lefering, R.; Neugebauer, E.; Beuth, J., König, D-P. Pulverer, G. Central Venous Catheters and Bloodstream Infection. JAMA. 2000. 28, 477.

Comparing zones of inhibition against *pseudomonas aeruginosa* that was cultured from Ruth A. Morrison (# 421222) as imparted by (a) the cook triple lumen polyurethane cvc, and (b) Gendine-impregnated triple lumen polyurethane.

Origin	Cook	Gendine
Peripheral	4 mm	13 mm
CVC	4mm	13mm

Brilliant Green

BG in DCM & MeOH

	MRS	MRSA ₂₀₆₆		PS ₄₂₀₅		C. Parap. ₁₋₁₀₀₋₀₀₂₂	
	DCM	MeOH	DCM	MeOH	DCM	MeOH	
PVC	30:30	16:17	0:0	0:0	27:27	11:11	
Si	0:0	0:5	0:0	0:0	0:0	0:0	
PU	20:20	18:18	0:0	0:0	25:25	21:21	
Silk	8:9	10:11	0:0	0:0	0:0	7:7	

BG+ T in DCM and Acetone

	MRS	MRSA ₂₀₆₆		PS ₄₂₀₅		C. Parap. 1-100-0022	
	DCM	Me ₂ CO	DCM	Me ₂ CO	DCM	Me ₂ CO	
PVC	18:18	18:19	0:0	0:0	13:13	11:11	
Si.	16:23	17:18	0:0	0:0	13:13	12:12	
PU	20:20	24:24	0:0	0:0	17:17	12:12	
Silk	12:12	10:11	0:0	0:0	8:8	7:7	

BG+CHX-/DCM

	MRS	MRSA ₂₀₆₆		PS ₄₂₀₅		0-1-100-0022
	DCM	MeOH	DCM	MeOH	DCM	MeOH
PVC	23:23	21:21	18:18	0:0	18:18	21:22
Si	14:17	15:14	8:9	0:0	18:18	6:8
PU	19:21	16:16	14:14	9:9	17:17	15:15
Silk	11:11	12:13	4:4	2:2	6:7	10:10

slightly soluble. Immersed for 24 h.

Brilliant Green (BG) with & without Chlorhexidine (CHX)*

	M	MRSA ₂₀₆₆		PS ₄₂₀₅		rap. ₁₋₁₀₀₋₀₀₂₂
	BG	BG*·CHX	BG	BG+CHX	BG	BG+CHX
PVC	30:30	23:23	.0:0	18:18	27:27	18:18
Si	0:0	14:17	0:0	8:9	0:0	18:18
PU	20:20	19:21	0:0	14:14	25:25	17:17
Silk	8:9	11:11	0:0	4:4	0:0	6:7

^{*}From DCM

Brilliant Green + Gentian violet, 1:1

	MRS	MRSA ₂₀₆₆		PS ₄₂₀₅		1-100-0022
	DCM	MeOH	DCM	MeOH	DCM	MeOH
PVC	28:28	17:17	13:13*	0:0	29:29	13:13
Si	8:9	0:0	0:0	0:0	0:0	0:0
PU	24:24	20:20	11:12*	0:0	26:26	17:17
Silk	7:7	7:7	0:0	0:0	0:0	0:6

^{*}Bacteristatic

Durability of Polyurethane & Silicone Impregnated with GV·CHX

Zones of Inhibition against MRSA2066 after x Days Incubation in Human Serum

	Catheter	Day=0	Day=3	Day=7	Day=10	Day=17	Day=31	Day=46	Day=61
L	PU	21:21	18:18	16:16	15:14	14:14			
[Si	21:21	18:18	15:15	11:11	11:11			


Durability of GV-CHX-coated Silicone UT Catheter

Zones of Inhibition against EN_{3836(VRE} and E. Coli ₃₂₂₆ after x Days Incubation in Human Urine

Organism	Day=0	Day=3	Day=7	Day=14	Day=21	Day=35	Day=42	Day=49
EN ₃₈₃₆	23:23	20	19					
E. Coli	18:18	15	13					

Day₀ against MRSA₂₀₆₆ = 26:27; against Ps = 18:18; against C. parap. = 24:25

	0	3	7	10	14	17	21	28
PU	21	18	16	14.5		14		
Si -	21	18	15	11		11		
ET-PVC	28				23	•	22	22.5



 30
 35
 43
 45
 50
 57
 61
 63


 11
 8
 6.5
 6.5

 7.5
 0
 0
 0

 22
 20.5
 22
 22
 19
 .21

77 5.5

0 7 10 14 21 28 35 43 MRSA 28 23 22 22.5 22 20.5



 50
 57
 63
 70
 77

 22
 22
 19
 21
 20

Durability of Gendine-Impregnated Devices

I. Durability of Polyurethane & Silicone Impregnated with GV-CHX

Zones of Inhibition against MRSA2066 after Incubation in Human Serum.

LOUICS OF	LIMIDICION	against II.	TACO 27 2000 0	arter aneu	Dation in	TIMMEN C	CI WAIA.	
Catheter	Day=0	Day=3	Day=7	Day=10	Day=17	Day=30	Day=45	Day=61
PU	21:21	18:18	16:16	15:14	14:14	11:11	8:8	6:7
Si	21:21	18:18	15:15	11:11	11:11	7:8	0:0	0:0
	Day=75							
PU	5:6							
Si	ND							

ND = Not done

II. Durability of Gendine-Coated Silicone UT Catheter

Zones of Inhibition against EN3836(VRE) and E. Coli 3226 after Incubation in Human Urine.

Organism	Day=0	Day=3	Day=7	Day=14	Day=21	Day=28	Day=35	Day=42	Day=49
EN3836	23:23	20	19:19	17:15	13:15	14:13	13:15	15:10	11:17
E. Coli	18:18	15 .	13:13	13:12	12:12	12:11	11:11	11:11	11:11
	Day=56	Day=63	Day=70	Day=78	Day=85.	Day=92	Day=99	Day=106	
EN ₃₈₃₆	Day=56 14:17	Day=63 14:18	Day=70 14:19	Day=78	Day=85	Day=92 17:17	Day=99 15:15	Day=106	

Day=0 against MRSA₂₀₆₆ = 26:27; against Ps₄₂₀₅ = 18:18; against C. parap.₁₋₁₀₀₋₀₀₂₂ = 24:25

III. Durability of Gendine-Coated ET PVC Tube

Zones of Inhibition after Incubation in Human BAL.

Organism '	Day=0	Day=7	Day=10	Day=14	Day=21	Day=28	Day=35	Day=43	Day=50
PS ₄₂₀₅	20:20	11:11	0:0						
MRSA ₂₀₆₆	28:28	ND	ND	23:23	22:22	22:23	22:22	20:21	22:22
	Day=57	Day=63	Day=70	Day=77					
MRSA ₂₀₆₆	22:22	19:19	21:21	20:20					

ND= Not done

IV. Duarability of Gendine-coated Polyurethane against $MRSA_{2066}$ using methanol as impregnating solvent

Day=0	Day=3	Day=7	Day=10	Day =14	Day=35	Day=46	Day =60	Day=91	
16:16	13:13	ND	12:13	14:14	11:11	10:11	9:9	7:6	




Figure 1. Comparison in durability between DCM and Methanol methods for impregnating polyurethane with Gendine upon incubation at 37° in human serum

Durability of polyurethane impregnated with GV (Gentian Violet).

About 30 one-centimeter sterilized double-lumen catheters (beige) were immersed in a solution of 2 g of GV in 60 ml MeOH (methanol) for 2 h. Catheters were removed from the solution, and traces of solution were removed from the lumen, then allowed to dry over night under the hood. The catheters were washed with distilled water (by shaking the catheters with water in a tube) till the washings were colorless or faint after which the catheters were allowed to dry for at least 4 h. The impregnated catheters were placed in a tube and covered with human serum (Sigma, # S7023), and allowed to stand for the appropriate period @ 37.5° C. Serum was replaced each time catheters were removed. The serum-soaked impregnated catheters were allowed to dry for at least 4 h, then embedded in agar plates (MH II) streaked with appropriate microorganism, and the plates were incubated for 24 h (corrected after 48 h). The resulting zones of inhibition are given below along with the control using cook catheters (impregnated with minocycline-rifampin).

Durability of polyUrethane Impregnated with GV (pUGV).

Zones of Inhibition against MRSA...

Zoncs or	mmonue	on agam	IST MILES	™2066∙				
Catheter	Day=0	Day=3	Day=7	Day=10	Day=17	Day=31	Day=46	Day=61
pUGV	24:24	20:20	19:19	21:21	16:16	15:13	12:12	12:10
Cook	22:22	17:16	20:20	21:21	29:30	11:11	7:7	0:0
-	Day=75		Van den en it in Nobel den in	-	and the second second	ATTRIA COARTERIOR	ACCUMENTATION OF THE PARTY OF T	S-100 C-200 C-
1	Day-10							
pUGV	9:10							

^{*}Cook catheter as control

Dr. Raad...

- 1. Results of doubling the concentration are given in Table 4.
- Table 5 summarizes results for GV+ PCMX using DCM. Notice the zones for silicone.
- The SIC catheters from Jim Yardly arrived yesterday. I have not done the experiments. I am trying to get hold of Al because I need the volumeter.
- Obtained similar zone of inhibition against MRSA for the arrow using both commercially available Mueller Hilton and TSA agar plates, where the zone are 18 mm.

Table 4. Comparison Between Zones at Concentrations 1x and 2x

	MRSA ₂₀₆₆		PS	3681	C. Parap.1-100-0022		
Catheter	1x	2x ·	1x	2x	1x	2x	
ET-PVC	24:26	28:29	15:16	20:21	25:29	28:29	
SI	12:13	13:13	0:0	0:0	14:15	13:13	
PU	20:20	21:21	8:10	12:12	20:21	25:22	
Silk	15:15	17:18	0:0	0:0	12:14	18:22	

¹x refers to 2.25 mmol of GV PCMX dissolved in 30 ml MeOH. 2x refers to 4.50 mmol.

Table 5.

	GV* PCMX' in DCM						
Catheter	MRSA	Ps	C. Parap.				
ET-PVC	25	19	25				
Si	20	12	18				
PU	19	15	21				
Silk	16	0:0	13				

^{*}Catheters are dipped in the solution until swelling becomes visible (few minutes).

Best wishes.

Nabeel

Cholestin 2 antippendannal abx.
polymixin }
Chloramine
Chlorhexidine - thymal varnish (Cervilec
acetiz acid
Zinc Chloride

Na hypochlorite

methyl iso thiazolone

India Gentian Violet improg silicone

DV Raed

LTI, Whil. assoc.

Condo. Gentian Violet

Frohim Fredical Research

get article

works less; ethane, methane

works best methylene Chlor. de

100

* Danna 212 733 6632 ned for #

Violet improg silicone


Works less i ethane, methane


Works best methylene Chlor. de

100

* Danna 212 733 6632 ned for #

	Icthammal (used: glycerol in ear chops)
, , (6)	stoph aureur ZI = 18 mm
	Strep pyagenes ZI = 23 mm
	nos gactivity against proteur mirabilis,
on Land	inhibited)
60°x	B. arrugineta, E. Cali (C. albitans, weakly inhibited) active against (T. Darg inhibited) hes anti-inflamm action
	I
	Mercurochrome (and in fungal ear inject. Ras antifungal effect: Aspergillus niger
	Ras antifungal effect: "Aspergillus niger
	A Flavas
	A. fumi galis
Committee of the second	Candida
	A. flavas A. fumigatus Candida Mucor
(3)	10-05/ Standard glacont - 0.025% benzal Konium
	Chloride + 4 / benzy/ alcohol
	VS. to / Poridone iodine
y	VS. 10 / Povidone iodine 7 Betadine
	VS. to / Betadine - Betadine used as Adutions for Cuc, arter. Cath. Care
Y	VS. to providence iodine Betadine used as Solutions for cuc, arter. Cath. Care Chlerkexid soln. was supporior to povidence iodine
	VS. to / Bridone iodine 7 Betadine used as Abutions for CVC, arter. Cath. Care 6 Chlerkxid-Solan was supporior to povidone iodine Bragainst G bacteria
	VS. to providence iodine Betadine used as Adutions for CVC, arter. Cath. Care Chlerkexid-Soln. was surperior to povidence iodine against Gr (1) bacteria Chlorkex seln. was nonsignificantly superior agains
	VS. to / Bridone iodine 7 Betadine used as Abutions for CVC, arter. Cath. Care 6 Chlerkxid-Solan was supporior to povidone iodine Bragainst G bacteria





may have an anti-inflammator ction through inhibition of ction through in activated DNA synthesis in activated

(6) (B) Alkaloids palmatine Sangunarine

- Inhibit the multiplication of bacteria, fung. and viruses

Sanguinarine -> inhibits choline acetyltransferase

· Berberine, palmatine - active at the

 2 receptor Berberine & Sanguinine intercalate DNA inhibit DNA synthesis and reverse transcript Sanguinarine affects membrane permeability Berberine affect protein biosynthesis

13-hexylberberine

several 13-alkyl substituted analogs of berberine and prematine are active against staph aureur.

- 13-heryl perberine 2 8X as active as 13-heryl palmatine Kanamycine sulfate

(12) Hovidone - iodine iodine and iodophors efficacions against meth reasist-Staph aureus (MRSA) Enterococcus, i no develop. J resistance excellent local tolerability of betaisodona preparations

		•			
		Povidone-iodine Killed X-malt	, Na hypo ophiha . S	chlorite marcescen	~ ~
		* Chlorhexidune	0.2%		
	•	didn't Kill . X. malt		narcescens (Capter 10 min
	. ,	* BenzalKoniun			
		1 Killed X m	altophilia		
	-	(0.1% was nee	ded to Kill	Simarces	Cers)
	ene hammy	* Tego = 51	(3)		
		0.02% Killed			
	(E)	Peroxyacetic	acid	(local a	ntiseptic)
		•			
		•		The second secon	
•	**	rg r	of the comment	name has sens " this see absolute	

Suture

Initial expriments: PCMX in MeOH + 50% aq. NaOH added to GV in MeOH, stirred for 1 h. Then evaporated solvent. Dissolved residue in DCM

Catheter	MRSA	Ps	C. Parap.
ET-PVC	23:22:25	12:0: <u>19</u>	22:25:25
Si	16:16:20	12:5:12	17:17:18
PU	20:21:19	16:16: <u>15</u>	23:25:21
Silk	12:12:16	0:0:0:0	12:12:13

¥

GV ⁺ ·⁻OH, prepa	ared by adding 50% aq. evaporating solvent, diss	NaOH to GV in MeO solving residue in DC	H, stirring for 1 h,
	MRSA	Ps	C. Parap.
PVC	20:20	0:0	17:17
Si	16:16	5:7	14:14
PU	20:20	14:15	25:25
Suture	13:14	0:0	12:12

GV+ OH, prepared by adding 50% aq. NaOH to GV in water, stirring for 1 h, evaporating solvent, dissolving residue in DCM. C. Parap. MRSA Ps PVC 23:23 0:0 1 22:22 Si 9:9 18:19 13:16 PII 23:23 17:17 21:21 Suture 17:17 3:4 13:14

GV+ OH PCMX, prepared by adding 50% aq. NaOH to GV in water, stirring for 1 h, evaporating solvent, dissolving residue in DCM. Then added PCMX MRSA Ps C. Parap. PVC. 22:22 0:0 20:20 Si 21:20 6:10 17:18 PU 25:25 14:14 24:26 Suture 15:16 3:3 13:14

GV+-PCMX, prepared by adding methanolic sodium methoxide to PCMX in MeOH. then the resulting mixuture is added to GV in MeOH, stirring for 1 h, evaporating the solvent, then dissolving residue in DCM MRSA Ps C. Parap. **PVC** 21:20 0:0 20:20 Si 18:18 10:10 25:25 PU 21:21 15:15 29:29

0:0

12:13

3

13:13

GV+ OCH3, prepare	ed by adding methanolic evaporating the solvent,	sodium methoxide to	GV in MeOH, stirring	
101 1 11, 0		then dissolving residue	III DCIVI	
	MRSA Ps C. Para			
PVC	21:20	0:0	20:20	
Si	18:19	9:9	25:25	
PU	21:21	15:15	28:29	
Suture	15:15	0.0	10:13	

aprotion Who who a protion with och ?

Comparison between MeOH and DCM solutions of GV.

	GV in MeOH			GV in DCM		
Catheter	MRSA	Ps	C. Parap	MRSA	Ps.	C. Paráp
ET-PVC	20-21	0-0	18-19	31-27	18-18	28-31
Si	7-8	0-0	0-0	10-10	0-0 🎷	0-7
PU	32-32	24-26	31-32	24-24	18-18	25-25
Silk	10-11	0-0	0-0	9-9	0-0	0-0

24

_

7.

1/27/00

Dr. Raad..

This report is to update you on where we are with Gentian violet & PCMX.

I have not received Clofoctol. I'll try to recover what we had from the first trials.

Nabeel

1/27/00

Summary of work with GV PCMX

Recall that PCMX Na+ was added to GV in methanol, and residue resulting from evaporation of methanol was dissolved in DCM. Table 1 shows zones of inhibitions obtained from this first attempt with DCM.

	Table 1.			
		GV [*] ·PCMX [*] in DCM		
. 1.41	Catheter	MRSA	Ps	C. Parap.
Vinyl Chlor	ET-PVC	25	19	25
licone	Si	20	12	18
wrethane	PU	19	15	21
suture	Silk	. 16	0:0	13

aichlorometham methylen chlorich Mp-Chloro-methyl3,5-divnethyl-Cxylenol

The experiment was repeated as follows:

Sodium hydroxide, 0.59 ml of 50% NaOH, was added to 1.15 g (7.35 mmol) of PCMX in 35 ml MeOH. The resulting solution was added dropwise to a solution of GV (3 g; 7.35 mmol) in 150 ml MeOH, and the resulting solution was stirred at ambient conditions for 1 h. The precipitate was filtered under vacuo. The filtrate was placed under the hood over night, allowing the solvent to evaporate. The resulting residue (3.279 g) was used without purification, of which 1.1 g was dissolved in 30 ml DCM for impregnating catheters. Results are given in Table 2.

Table 2.				
	GV ⁺ ·PCMX [*] in DCM			
Catheter	MRSA	Ps	C. Parap.	
ET-PVC	22:23	12:0/	22:25	
Si	16:16	12:5	17:17	
PU	20:21	16:16	23:25	
Citt	42.42	0.0	40.40	

PU was immersed for 20 sec. Silk suture was immersed for 2 h.

In order to optimize the procedure, it is important to examine the chemistry closely. Careful examination of GV reveales that it is not a true quaternary amine. Quaternary amines have a net positive charge localized on the nitrogen because their lone-pair electrons are involved in covalent bonding with a nucleophilic center. An example is tetramethylammonium iodide (Equation 1).

Quaternary ammonium compounds (Quats) are used as phase-transfer catalysts (PTC) for solvating organic salts in organic media by increasing the lipophilicity of the organic salt via formation of an ionic compound with the PTC (Figure 1).

The case is different for Gentian violet. It is not a true quat. It is an amine with a benzylic carbonium center, which is stabilized by resonance through, in addition to the phenyl rings, delocalizing of the lone-pair of the nitrogens via resonance.

As a result, the positive charge is not localized on the nitrogen. Briefly, and depending on the resonance hybride of GV, it can form ionic compounds with organic salts, and depending on the basicity of the nucleophile, can form covalent adducts at the carbonium center. In other words, the presence of polar protic solvents can compete with the nucleophile.

In the case of PCMX, water, methanol can react with GV at the carbonium center in the presence of a proton acceptor. The phenoxide derivative of PCMX is not a strong nucleophile, but can accept a proton from a hydronium ion. This implies multiple product formation when PCMX is added to a solution of GV in methanol and the presence of water. In other words, the following can result from addition of PCMX to GV in alcoholic aqueous media:

- 1) GV+PCMX.
- ⊤2) GV⁺·OH.
 - 3) GV⁺· OCH₃⁻. ~>

The sure way to shed more light on the results on hand, preparing each of these possible reagents and testing their impregnating ability and zones of inhibitions will help shed more light onto the results on hand.

The following results summarize work to date.

Table 3.

	GV⁺∙⁺OI	I/MeOH	
	MRSA	Ps	C. Parap.
PVC	20:20	0:0	17:17
Si	16:16	5:7	14:14
PU	20:20	14:15	25:25
Suture	13:14	0:0	12:12

Table 4

	GV ⁺ · OC	H ₃ /MeOH		
MRSA Ps C.				
PVC	21:20	0:0	20:20	
Si	18:19	9:9	25:25	
PU	21:21	15:15	28:29	
Suture	15:15	. 0:0 .	10:13	

Table 5

ore o.				
	GV ⁺ ·PCMX/MeOH			
	MRSA	.Ps	C. Parap.	
PVC	21:22	-0:0	20:20	
Si	18:18	10:10	. 25:25	
PU	21:21	15:15	29:29	
Suture	12:13	0:0	13:13	

Two more experiments in this series are undergoing.

- -The experiment employing acetonitrile as a solvent was accidentally spilled.
- -I'll set up the rifampin experiment as soon as possible.